

# 全球血糖機參考機種

香港 (Hong Kong)

## 重要資訊:

這個目錄是參考網際網路，文獻，百特內部的研究，政府機關，試紙宣傳單，安全性報告，和來自產品製造商的訊息。我們已經盡力提供準確和最新的資訊。由於百特不生產這些血糖機或試紙片，所以不能保證這些資訊的準確度。當EXTRANEAL (Icodextrin) 腹膜透析液與血糖機共同使用前，請與製造商聯繫您的血糖機或試紙的兼容性資訊。

- EXTRANEAL (Icodextrin) 腹膜透析液包含icodextrin的成分。麥芽糖是icodextrin 的代謝產物之一，麥芽糖會干擾這些血糖機或血糖試紙，如使用某些特定的血糖機或血糖試紙，可能會有錯誤的血糖值被誤判的結果。
- 錯誤的血糖判讀，可能掩蓋實際上是血糖過低的情況，而導致誤判為血糖過高。使用EXTRANEAL (Icodextrin)腹膜透析液的病人使用上述血糖機或試紙，血糖判讀為正常值時，可能掩蓋了低血糖的情形。此將導致病人或醫護人員沒有採取適當的措施使血糖值回復正常範圍，誤判病人為高血糖，則可能讓病人接受過多的胰島素治療。這兩種情況都有可能威脅到病人安全，包括喪失意識，昏迷，神經損傷或死亡。
- 透析中心及病人都應該直接與您血糖機或試紙製造商聯繫，其測試方式是否會被EXTRANEAL (Icodextrin)腹膜透析液或麥芽糖干擾血糖值，或是只會偵測到葡萄糖。請詳閱血糖機或試紙內附的產品說明。
  - 以下資訊僅供參考。並不意指您該選擇何者廠牌血糖機或試紙。
  - 測試方式請參考以下表格。**使用EXTRANEAL (Icodextrin)腹膜透析液的病人只能使用具有葡萄糖專一性的血糖機或試紙來測量血糖。**各試紙與相配的血糖機也列舉在以下廠商名單中。如果您是使用其他試紙，請與製造商聯繫來確認您血糖試紙是否具有葡萄糖專一性。以下資訊為血糖機廠商主要及較常見的血糖機或試紙。
  - 此並非完整的血糖機及試紙製造商/供應商名單(更新日期為二零一零年六月)。若您使用的血糖機或試紙不在此，不代表其不可使用於EXTRANEAL (Icodextrin)腹膜透析液的病患。如有任何疑問都請與製造商聯繫以得到更進一步的資訊。百特沒有責任定期更新此名單。

## 血糖機及試紙製造商/供應商參考

更新日期為二零一一年六月

血糖機廠牌/商品名	與Extraneal (Icodextrin) 兼容性 (葡萄糖專一性)	測試方法*	製造商/供應商
Boots Omron HEA-214 Optium Optium Easy OptiumEZ Optium Xceed Optium Xido Precision PCx Precision QID <sup>B</sup> Precision Xceed Precision Xceed Pro Precision Xtra Precision Xtra OK ReliOn Ultima Sof-Sense Sof-Tact <sup>C</sup> TrueSense	是 是 是 是 是 是 是 是 是 <sup>1</sup> 是 是 是 是 是 是 是 是 是 是 是	GDH-NAD GDH-NAD GDH-NAD GDH-NAD GDH-NAD GDH-NAD GDH-NAD GDH-NAD, GO <sup>1</sup> GO GDH-NAD GDH-NAD GDH-NAD GDH-NAD, GO <sup>1</sup> GDH-NAD GDH-NAD GDH-NAD	Abbott Diabetes Care <a href="http://www.abbottdiabetescare.com">www.abbottdiabetescare.com</a> (852) 2806 4488





參考資料:

- <sup>A</sup> Baxter study 32386 (McGaw Park) Determination of potential interference of icodextrin and its metabolites on human blood glucose measurement using Accu-Chek compact and Advantage systems.
- <sup>B</sup> Baxter report REP-NIV-RE-366, (Nivelles) Evaluation of potential interference in blood glucose determination (measured with enzymatic methods) for patients treated with icodextrin.
- <sup>C</sup> Baxter report RD-01-RE-233 (McGaw Park) Evaluation of potential interference in the blood glucose test kit MediSense Softact with icodextrin and its metabolites.
- <sup>D</sup> Baxter report 33541 (McGaw Park) Determination of potential interference of icodextrin and its metabolites on human blood glucose measurement using chosen glucometer-Glucocard X-Meter (Arkray)
- <sup>E</sup> Baxter report Interim 1, 33541 (McGaw Park) Determination of potential interference of icodextrin and its metabolites on human blood glucose measurement using chosen glucometers.

請參考EXTRANEAL (icodextrin) 腹膜透析液的完整處方說明。  
Baxter and EXTRANEAL™ 為百特醫療用品有限公司之註冊商標。  
其他所標示商標及產品均歸屬其製造公司