

WARNING

Potential for Incorrect Blood Glucose Reading

Is using **EXTRANEAL (Icodextrin 7.5%)**
peritoneal dialysis solution

0508

ONLY use glucose monitors and test strips that are glucose-specific. These methods are common in clinical laboratories. Contact the manufacturer of the glucose monitors and test strips to determine the method that is used.

Do not use glucose monitors or test strips that utilize the enzyme glucose dehydrogenase pyrroloquinolinequinone (GDH-PQQ) or glucose-dye-oxidoreductase. This may result in a falsely elevated blood glucose reading in patients using **EXTRANEAL (Icodextrin 7.5%)**. This may mask true hypoglycemia or lead to the erroneous diagnosis of hyperglycemia, leading to life-threatening events.

For further information, visit www.glucosafety.com

Baxter BAXTER and EXTRANEAL are trademarks of Baxter International Inc.

102.797mm x 40mm