

WARNING

Potential for Incorrect Blood Glucose Reading

Patient Name

Is using EXTRANEAL (icodextrin)
peritoneal dialysis solution

Baxter

BAXTER and EXTRANEAL are trademarks
of Baxter International Inc. AL06311A 06/08

ONLY use glucose monitors and test strips that are glucose-specific. These methods are common in clinical laboratories. Contact the manufacturer of the glucose monitors and test strips to determine the method that is used. DO NOT use glucose monitors or test strips that utilize the enzyme glucose dehydrogenase pyrroloquinolinequinone (GDH PQQ) or glucose-dye-oxidoreductase. This may result in a falsely elevated blood glucose reading in patients using EXTRANEAL (icodextrin). This may mask true hypoglycemia or lead to the erroneous diagnosis of hyperglycemia, leading to life-threatening events.

For further information call 1-888-736-2543 option 1.

For full prescribing information on EXTRANEAL and additional glucose monitor information visit www.glucosesafety.com.